

## CryoBanker® Cell Freezing Medium, Serum Free

Cat. No. CB-50, CB-100

Storage Temperature 4°C

CRYOBANKER is a serum-free cell cryopreservation medium with unique complete formulation for broad spectrum of mammalian cell cultures. The unique formulation allows for stable cryopreservation and high viability after freeze-thaw procedure, without potential risk of contamination by serum. Optimum for serum-free cultured cells and peptide/protein expressing cells cryopreservation.

### 1. Characteristics

- High cell viability,
- No risk of contamination by serum derived pathogens
- Batch-to-batch stability.
- Complete formulation,
- No further preparation,
- Direct freezing at -80°C,
- No programmable freezer required, etc.

### 2. Contents

- Serum-free formulation.
- Contains DMSO.

Cat. No.	Quantity
CB-50	50ml
CB-100	100ml

### 3. Storage/Stability Conditions

The medium is shipped on cool packs. Upon receipt, store at 2~8°C. The medium is stable until the expiration date stated on the label.

### 4. Protocols

#### 1. Freezing

For optimum results, cells for cryopreservation should be in log phase of growth. Similar or standard freezing protocols may be substituted.

- 1) Examine and make sure the cell culture free of contamination, in healthy situation and proper confluency, etc.
- 2) Perform a cell count to determine the viability of cells.
- 3) Gently pellet the cells by centrifugation (3 - 5 minutes at 1,000~2,000rpm, 4°C). Remove the supernatant by using an aspirator.
- 4) Gently suspend the cells with CRYOBANKER® cryopreservation medium (1 ml for  $5 \times 10^5$ - $6$  cells).
- 5) Dispense the cell suspension in 1ml aliquots to cryopreservation vials that have been labeled with the cell line name, cell concentration, passage date and other essential information.
- 6) Place the vials directly -80°C for storage. If necessary, transfer the frozen vials to a liquid nitrogen storage tank after the vials have been frozen for at least 24 hours.

#### II. Thawing

- 1) Remove the frozen cell from storage and quickly thaw in a 37°C shaking water bath.
- 2) Immediately dilute and gently mix each 1ml of cells with 10ml of complete cell culture medium.
- 3) Gently pellet the cells by centrifugation (3-5 minutes at 1,000 ~ 2,000rpm, 4°C). Remove the supernatant by aspirator.
- 4) Gently suspend the cells with appropriate volume of complete cell culture medium. And plate in a culture flask.
- 5) Continue the further culture procedures according to standard protocols.

### 5. Quality Control

Quality Control TESTS	SPECIFICATION
Endotoxin tested	Not more than 6 EU/ml
Mycoplasma tested	Not detected
Cell culture tested	Hela cell
Sterile-filtered	0.2um

### 6. Applicable cell lists

A10	CHO	ELM-D	Kato-III	NIH3T3	Vero
A2781	CHO-K1	HeLa cell	KM12-LX	P388	WEHI3B
Ac2F	Colo203	Hep3B	L5178Y	PC12	WiDr
AfT20	COS1	Hepal-6	L929	Raji	293
Ba/F3	COS7	HepG2	LM	RAW264.7	293T
BHK-21	CTL-2	HL-60	LNCap	sf9	32D
C2C12	DLP-1	HT-2	MCF-7	SK-N-MC	3LL
C3H10T1/2	DT40	Huh-7	MDCK	SN12C	
Caco-2	DU145	Jurkat	Molt-4	THle3	
CHL/1U	EJ-1	K-562	NCI-H441	T24	

### 7. Related Products

*BioMycoX*® Mycoplasma Detection Kit  
Cat. No. D-25, D-50, D-100

*BioMycoX*® Mycoplasma Elimination Kit  
Cat. No. E-01 (1 application)

*BioMycoX*® Mycoplasma Prevention Spray  
Cat. No. P-1000 (1 L)